Oxygen Catheter

DESCRIPTION

The Oxygen Catheter is a sterile, single-use tube. It is used to deliver oxygen to a patient. The Oxygen Catheter is made from PVC in medical grade. The Oxygen Catheter have two models, with fixed compress and plain type. The fixed compress type consist of connector, tubing and sponge distal end. The plain type consist of connector and tubing.

This product can be used on patient who need low concentration oxygen supply.

FEATURES

- Suitable for easy and efficient administration of oxygen
- Soft and smooth open distal end facilitates non-traumatic insertion
- Funnel connector at proximal end

INTENDED USE

The Oxygen Catheter is intended for direct administration of oxygen via nasopharyngeal route.

CLINICAL BENEFITS

The Oxygen catheter can deliver the oxygen to patients via nasopharyngeal route, the oxygen could enter into the body and maintain normal oxygen saturation.

The clinical benefit associated with the device is define as direct, it is intended used for Infant, child and Adult patients who need low flow of oxygen therapy.

PATIENT TARGET GROUP

Infant: 8Fr

Child: 8Fr, 10Fr

Adult: 12Fr, 14Fr and 16Fr

INTENDED USER

Must be operated by trained professionals.

INDICATIONS

If the following diseases are accompanied by hypoxemia without hypercapnia and oxygen flow is required to be 1-5 L/min, the Oxygen Catheter should be selected for oxygen therapy:

- Decreased vital capacity
- Cardiac insufficiency
- Dyspnea caused by various poisonings
- Coma
- Postoperative patients with massive hemorrhagic shock or craniocerebral diseases.

CONTRAINDICATIONS

- Complete obstruction of nasal meatus.
- Patients who are allergic to plastics are not recommended to use the Oxygen Catheter.

ADVERSE REACTIONS

- Possible headache or mucosa dryness.
- Possible skin irritation or rupture. Tight fixation is prohibited.

- Respiratory depression. Patients with type II respiratory failure should be given low concentration of oxygen at low flow (1~2 L/min) to maintain PaO₂ at 60mmHg.

DIRECTIONS FOR USE

- Check if the package is intact and make sure the product is within the expiry date.
- Select one nose and use wet cotton swab to clean the nose.
- Connect Oxygen Catheter to the oxygen supply system via an Oxygen Tubing.
- Adjust to suitable flow rate.
- Oxygen Catheter with fixed compress: carefully insert compress into one nostril.

Oxygen catheter without compress: carefully insert catheter via one nostril, down the nasopharyngeal route placing the catheter tip no further than the uvula. The insertion length is usually 2/3 of distance from the nasal tip to the earlobe. Secure catheter via a securement device.

PRECAUTIONS/ WARNINGS

PRECAUTIONS

- Choose the proper size before using it. Check the model and size on the product label to make sure the proper model and size is used.
- Adjust the oxygen flow rate first before oxygen supply. When stop oxygen supply, remove the oxygen tubing first and then switch off oxygen source.

WARNINGS

- Read all warnings and instructions before use. Improper use can result in serious or fatal illness or injury.
- This device can only be used by trained professionals.
- Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Since oxygen is a flammable substance, there should always be a “No Smoking” sign posted when a person is receiving oxygen.
- Keep the oxygen container away from open flames or heaters, such as candles, fireplaces, gas stoves, or hot water heaters.
- Be sure all connections are secure.
- Do not use if it's damaged or irregularly shaped.
- Do not use after the expiry date.
- For single patient use only.
- Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in injury, illness or death of the patient. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- Dispose of product and packaging in accordance with hospital administrative and/or local government policy.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the local competent authority of the user's place.

STORAGE CONDITIONS

- Store in a cool and dry place, keep away from sunlight.

- Protect product from moisture and excessive heat.
- Avoid prolonged exposure to ultraviolet, sunlight and fluorescent light.
- Store in manner preventing crushing.

The Key Performance Characteristics of Oxygen Catheter

- The Key performance: Resistance to flow

The pressure increment per meter length of pipe shall not exceed 2.0kPa when tested at the test flow rate corresponding to each size of product in the table below.

| Size(Fr) | Test flow(L/min) | Pressure increment per meter length of catheter(kPa) |
|----------|------------------|--|
| 8Fr | 1.0 ±0.3 | ≤200 |
| 10Fr | 2.0±0.3 | ≤200 |
| 12Fr | 3.5±0.3 | ≤200 |
| 14Fr | 4.0±0.3 | ≤200 |
| 16Fr | 4.0±0.3 | ≤200 |

- The Key performance: Bending resistance

Place the Oxygen Catheter in the environment temperature of 40 °C and 5 °C for 24h, repeated bending simple nasal oxygen tube 10 times, the catheter should not rupture, the inner wall of the lumen should not be bonded.

MEANING OF SYMBOLS ON PACKAGE

Date of manufacture

Use-by-date

Batch code

Do not re-use

Consult instructions for use

Do not use if package is damaged

Sterilized using ethylene oxide

Do not re-sterilize

Single sterile barrier system

Medical Device

Manufacturer

Authorized representative in the European Community

CE 0123

CE Marked
Product

Keep dry

Keep away from
sunlight

Importer

This way up

Catalogue
number

Fragile, handle
with care

Not contains or
presence of
phthalate.
(if applicable)

Temperature
limits:
-15°C~49°C

WELL LEAD MEDICAL CO., LTD.

Address: C-4 Jinhua Industrial Estate, Hualong, 511434 Panyu, Guangzhou, PEOPLE'S
REPUBLIC OF CHINA

E-mail: info@welllead.com.cn www.welllead.com.cn

Tel: +86 20 84758878

Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestrasse 80, 20537 Hamburg, GERMANY

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