

QSafe Laryngeal Mask



WELL LEAD MEDICAL CO., LTD.

Add: 47 Guomao Avenue South, Panyu Guangzhou, China

Marketing&Trading Office

Add: 19A Lancun Mansion, 1381 Dongfang Road, Shanghai 200127, P.R.China

Tel: +86-021-53897623 Fax: +86-021-53897615 Post Code: 200127

www.wellead.com

e-mail: info@wellead.com

wellead[®]

Intended Use

QSafe Laryngeal Mask is used for anesthesia, emergency, and resuscitation.

Features

Integrated bite block prevents
collapsion

Anatomically shaped tube
facilitates rapid placement

Soft anatomical cuff
provides good seal

Color-coding for
convenient size
identification

Drainage tube reduces
the risk of aspiration

Specification

Size	Patient Size	Max Cuff Inflation Volume	Max Gastric Tube	Cat. No.
3	Small adults 30-50 kg	30ml	14Fr	A05A023010
4	Medium adults 50-70 kg	40ml	14Fr	A05A024010
5	Large adults 70-100 kg	40ml	14Fr	A05A025010