


wellead

Enteral Nutrition

ENFit

The new international standard design for enteral feeding devices that complies with ISO 80369-3, which will provide better patient safety and avoid tubing misconnections. The new standard has changed the configuration of the feeding connectors and will not allow connectivity with the luer connectors or any other connectors for other clinical use. All Well Lead enteral feeding product series now are available with safer ENFit connectors to keep your patients properly connected all the time.




● PVC Stomach Tube

Features

- Made of medical grade PVC
- X-ray opaque line throughout the length
- With 4 side holes
- Available with open and closed distal end
- Available with funnel and ENFit connector

Size	Length (cm)	Cat. No. (With Open Tip & Funnel Connector)	Cat. No. (With Closed Tip & Funnel Connector)	Cat. No. (With Open Tip & ENFit Connector)	Cat. No. (With Closed Tip & ENFit Connector)
6Fr	120	N03A01061F	N03A01063F	N03A01062F	N03A01066F
8Fr	120	N03A01081F	N03A01083F	N03A01082F	N03A01086F
10Fr	120	N03A01101F	N03A01103F	N03A01102F	N03A01106F
12Fr	120	N03A01121F	N03A01123F	N03A01122F	N03A01126F
14Fr	120	N03A01141F	N03A01143F	N03A01142F	N03A01146F
16Fr	120	N03A01161F	N03A01163F	N03A01162F	N03A01166F
18Fr	120	N03A01181F	N03A01183F	N03A01182F	N03A01186F
20Fr	120	N03A01201F	N03A01203F	N03A01202F	N03A01206F
22Fr	120	N03A01221F	N03A01223F	N03A01222F	N03A01226F



● Silicone Stomach Tube

Features

- Made of 100% medical grade silicone
- X-ray opaque line throughout the length
- With 4 side holes
- Markings from the distal end every 5 cm
- Available with funnel and ENFit connector

Size	Length (cm)	Cat. No. (With Funnel Connector)	Cat. No. (With ENFit Connector)
8Fr	120	N03A020810	N03A020860
10Fr	120	N03A021010	N03A021060
12Fr	120	N03A021210	N03A021260
14Fr	120	N03A021410	N03A021460
16Fr	120	N03A021610	N03A021660
18Fr	120	N03A021810	N03A021860
20Fr	120	N03A022010	N03A022060



PVC Feeding Tube

■ Features

- Made of medical grade PVC
- X-ray opaque line throughout the length
- Available with Luer and ENFit connector

Length (cm)	Size	Cat. No. (With Luer Connector)	Cat. No. (With ENFit Connector)
50	4 Fr	N03C01041F	N03C01042F
	5 Fr	N03C01051F	N03C01052F
	6 Fr	N03C01061F	N03C01062F
	8 Fr	N03C01081F	N03C01082F
	10 Fr	N03C01101F	N03C01102F
	12 Fr	N03C01121F	N03C01122F
	14 Fr	N03C01141F	N03C01142F
	16 Fr	N03C01161F	N03C01162F
	18 Fr	N03C01181F	N03C01182F
	20 Fr	N03C01201F	N03C01202F
80	4 Fr	N03C0104AF	N03C0104BF
	5 Fr	N03C0105AF	N03C0105BF
	6 Fr	N03C0106AF	N03C0106BF
	8 Fr	N03C0108AF	N03C0108BF
	10 Fr	N03C0110AF	N03C0110BF
	12 Fr	N03C0112AF	N03C0112BF
	14 Fr	N03C0114AF	N03C0114BF
	16 Fr	N03C0116AF	N03C0116BF
	18 Fr	N03C0118AF	N03C0118BF
	20 Fr	N03C0120AF	N03C0120BF
120	4 Fr	N03C0104CF	N03C0104DF
	5 Fr	N03C0105CF	N03C0105DF
	6 Fr	N03C0106CF	N03C0106DF
	8 Fr	N03C0108CF	N03C0108DF
	10 Fr	N03C0110CF	N03C0110DF
	12 Fr	N03C0112CF	N03C0112DF
	14 Fr	N03C0114CF	N03C0114DF
	16 Fr	N03C0116CF	N03C0116DF
	18 Fr	N03C0118CF	N03C0118DF
	20 Fr	N03C0120CF	N03C0120DF



ENFit



Silicone Feeding Tube

■ Features

- Suitable for gastric feeding and enteric feeding
- Made of medical grade silicone material
- The ENFit connector engineered to reduce risk of tubing misconnections.
- The Y-port ENFit connector enabling simultaneous uninterrupted tube feeding and medication administration.
- X-ray opaque line throughout the length
- Imprinted with FR size, length, and depth scale for quick and precise placement
- Equipped with 4 side holes ensure continuous flow and minimize occlusion
- Transparent tube material for visual inspection for administration of nutrition and liquid



Size	Length (cm)	Cat. No. (With Guide Wire)	Cat. No.
4Fr	40	N03A120440	N03A110440
4Fr	75	N03C120475	N03C110475
4Fr	125	N03C120412	N03C110412
6Fr	55	N03A120655	N03A110655
6Fr	75	N03A120675	N03A110675
6Fr	95	N03A120695	N03A110695
6Fr	125	N03A120612	N03A110612
8Fr	55	N03A120855	N03A110855
8Fr	75	N03A120875	N03A110875
8Fr	95	N03A120895	N03A110895
8Fr	125	N03A120812	N03A110812
8Fr	140	N03A120814	N03A110814
10Fr	95	N03A121095	N03A111095
10Fr	125	N03A121012	N03A111012
10Fr	140	N03A121014	N03A111014
12Fr	95	N03A121295	N03A111295
12Fr	125	N03A121212	N03A111212
12Fr	140	N03A121214	N03A111214
14Fr	95	N03A121495	N03A111495
14Fr	125	N03A121412	N03A111412
14Fr	140	N03A121414	N03A111414
16Fr	95	N03A121695	N03A111695
18Fr	95	N03A121895	N03A111895
20Fr	95	N03A122095	N03A112095

◆ PU feeding Tube

■ Features

- Suitable for gastric feeding and enteric feeding
- Made of medical grade Polyurethane (PU) material
- The ENFit connector engineered to reduce risk of tubing misconnections.
- The Y-port ENFit connector enabling simultaneous uninterrupted tube feeding and medication administration.
- Available with Weighted and Standard (Bullet) tip
- X-ray opaque line throughout the length
- Imprinted with FR size, length, and depth scale for quick and precise placement
- Multiple side holes ensure continuous flow and minimize occlusion
- Transparent tube material for visual inspection for administration of nutrition and liquid



Bullet tip



Weighted Tip



● Bullet Tip

Size	Length (cm)	Cat. No. (With Guide Wire)	Cat. No.
4Fr	40	N03A080440	N03A070440
4Fr	75	N03A080475	N03A070475
4Fr	125	N03A080412	N03A070412
6Fr	55	N03A080655	N03A070655
6Fr	75	N03A080675	N03A070675
6Fr	95	N03A080695	N03A070695
6Fr	125	N03A080612	N03A070612
8Fr	55	N03A080855	N03A070855
8Fr	75	N03A080875	N03A070875
8Fr	95	N03A080895	N03A070895
8Fr	125	N03A080812	N03A070812
8Fr	140	N03A080814	N03A070814
10Fr	95	N03A081095	N03A071095
10Fr	125	N03A081012	N03A071012
10Fr	140	N03A081014	N03A071014
12Fr	95	N03A081295	N03A071295
12Fr	125	N03A081212	N03A071212
12Fr	140	N03A081214	N03A071214
14Fr	95	N03A081495	N03A071495
14Fr	125	N03A081412	N03A071412
14Fr	140	N03A081414	N03A071414
16Fr	95	N03A081695	N03A071695
18Fr	95	N03A081895	N03A071895

● Weighted Tip


Size	Length (cm)	Cat. No. (With Guide Wire)	Cat. No.
4Fr	75	N03C100475	N03C090475
6Fr	95	N03C100695	N03C090695
8Fr	125	N03C100812	N03C090812
8Fr	140	N03C100814	N03C090814
10Fr	125	N03C101012	N03C091012
10Fr	140	N03C101014	N03C091014
12Fr	125	N03C101212	N03C091212
12Fr	140	N03C101214	N03C091214
14Fr	125	N03C101412	N03C091412
14Fr	140	N03C101414	N03C091414

Gastrostomy Feeding Tube

Gastrostomy feeding tube is designed for long term enteral nutrition. It is inserted through a small incision in the abdomen into the stomach. It is useful where there is difficulty with swallowing with patient.

■ Features

- Made of medical grade silicone material
- Gastric balloon for a secure and comfortable fix
- Rounded tip with open distal end
- Rounded skin disk can easily adjust to maintain proper tube position
- Graduated shaft allows for easy stoma depth measurement
- Available with Luer and ENFit connector
- Guide wire is available



Connector	Size	Cat. No. (With Guide Wire)	Cat. No.
Luer	12Fr/5mL	N03D011211	N03D011210
	14Fr/5mL	N03D011411	N03D011410
	16Fr/15mL	N03D011611	N03D011610
	18Fr/15mL	N03D011811	N03D011810
	20Fr/15mL	N03D012011	N03D012010
	22Fr/15mL	N03D012211	N03D012210
	24Fr/15mL	N03D012411	N03D012410
ENFit	12Fr/5mL	N03D011221	N03D011220
	14Fr/5mL	N03D011421	N03D011420
	16Fr/15mL	N03D011621	N03D011620
	18Fr/15mL	N03D011821	N03D011820
	20Fr/15mL	N03D012021	N03D012020
	22Fr/15mL	N03D012221	N03D012220
	24Fr/15mL	N03D012421	N03D012420



WELL LEAD MEDICAL CO., LTD.

Address: 47 Guomao Avenue South, Panyu, Guangzhou, China

Tel: +86-20-84758878 Fax: +86-20-84758224 Post Code: 511434

Marketing & Trading Office

Add: RM 805, No 55, 21 Long, Pusan Road, Shanghai, China

Tel: +86-21-58397623 Fax: +86-21-58397615 Post code: 200125

E-mail: info@welllead.com.cn

<http://www.welllead.com>

wl.Enteral Nutrition.v1.261008