


Ambulatory Infusion Pump



Additional PCA Bolus



Intermittent Mode



Timed Mode



ISO13485:2016  0123 NMPA Registered

WELL LEAD MEDICAL CO., LTD.
Tel: +86-20-84758878 E-mail: info@welllead.com.cn
www.welllead.com

wellead |  Medsure
make sure our life

EASY

Simple operation
Easy parameter setting

SAFE


Five-year service life
Complete alarm system

SMART

Real-time monitoring system
Multiple infusion modes
Smart delivery



Medsure 3000 series ambulatory infusion pump supports multiple innovative therapies.



Pain management



Long-term IV antibiotics



Chemotherapy



Home treatment

SPECIFICATIONS

Type:	3000F
Cassette sizes:	50mL,100mL,150mL,200mL,300mL
Power supply:	2 AA batteries; at least infusion 100 hours in 2mL/h
Occlusion level:	High,medium,low; range: 30kpa ~ 180kpa

DEHP Free

Infusion parameters:

Total volume	0.1-999mL
Loading dose	0.1-30mL
Continuous rate	0.1-120mL/h
Maximum dose	1-120mL/h
PCA dose	0.1-9.9mL
PCA lockout	1-120min
Intermittent time	1min-48hr59min
Intermittent dose	1-999mL
Dose duration	1s-59min59s
Timed rate	0.1-120mL/h
Time Settings	24 hours with minimum step 15minutes
Infusion accuracy	±6%
Service life	5 years

Complete alarm system:

Occlusion Alarm, Air Alarm, Cassette Detached, Low Battery, Infusion Nearly Completed, etc.

For details, please refer to the user manual.