

Dual Seal Laryngeal Mask with Gastric Access



WELL LEAD MEDICAL CO., LTD.

Address: 47 Guomao Avenue South, Panyu, Guangzhou, China

Tel: +86-20-84758878 Fax: +86-20-84758224 Post Code: 511434

Marketing & Trading Office

Add: RM 805, No 55, 21 Long, Pusan Road, Shanghai, China

Tel: +86-21-58397623 Fax: +86-21-58397615 Post code: 200125

E-mail: info@welllead.com.cn www.welllead.com

wellead[®]

Intended Use

The dual seal laryngeal mask with gastric access aims to minimize the risk of aspiration and provide high seal pressure with its dual seal design and gastric-access channel.

Features

- Gastric access channel for managing gastric content
- Stylet shapes anatomical curve, ensuring rapid and easy placement
- Medical grade silicone allows better biocompatibility
- Silicone cuff creates anatomical seal of the pharyngeal
- Disposable laryngeal mask reduces the risk of cross-contamination
- Pilot balloon shows mask size and provides tactile indication of inflation degree
- Bite block protects airway occlusion



Specification

Mask Size	Patient Size	Largest Size of Gastric Tube	Cat. No.
1.5	5-10kg	10Fr	A05B021510
2	10-20kg	10Fr	A05B022010
2.5	20-30kg	14Fr	A05B022510
3	30-50kg	16Fr	A05B023010
4	50-70kg	16Fr	A05B024010
5	70-100kg	18Fr	A05B025010