


Latex Foley Catheter Ordering Information

■ Features

- Made of latex
- Silicone coated
- Color coding for size identification
- Sterilized by EO
- For single use only
- Available with different types



wellead

WELL LEAD MEDICAL CO., LTD.

Address: 47 Guomao Avenue South, Panyu, Guangzhou, China
Tel: +86-20-84758878 Fax: +86-20-84758224 Post Code: 511434
Marketing & Trading Office
Add: RM 805, No 55, 21 Long, Pusan Road, Shanghai, China
Tel: +86-21-58397623 Fax: +86-21-58397615 Post code: 200125
E-mail: info@wellead.com.cn
<http://www.wellead.com>

Latex Foley Catheter Ordering Information

Size (Fr/Ch)	Color Code	Length (mm)	Balloon (ml)									
			3ml	3-5ml	5ml	5-10ml	10ml	5-15ml	15ml	15-30ml	30ml	30-50ml
2-way Pediatric												
6	Light Red	270	F01A060602	—	—	—	—	—	—	—	—	—
8	Black	270	F01A060802	F01A060803	F01A060804	—	—	—	—	—	—	—
10	Grey	270	F01A061002	F01A061003	F01A061004	—	—	—	—	—	—	—
2-way Standard												
12	White	400	—	—	F01A021204	F01A021205	F01A021207	F01A021206	F01A021208	F01A021209	F01A021210	—
14	Green	400	—	—	F01A021404	F01A021405	F01A021407	F01A021406	F01A021408	F01A021409	F01A021410	—
16	Orange	400	—	—	F01A021604	F01A021605	F01A021607	F01A021606	F01A021608	—	F01A021610	F01A021612
18	Red	400	—	—	F01A021804	F01A021805	F01A021807	F01A021806	F01A021808	—	F01A021810	F01A021812
20	Yellow	400	—	—	F01A022004	F01A022005	F01A022007	F01A022006	F01A022008	—	F01A022010	F01A022012
22	Violet	400	—	—	F01A022204	F01A022205	F01A022207	F01A022206	F01A022208	—	F01A022210	F01A022212
24	Blue	400	—	—	F01A022404	F01A022405	F01A022407	F01A022406	F01A022408	—	F01A022410	F01A022412
26	Pink	400	—	—	F01A022604	F01A022605	F01A022607	F01A022606	F01A022608	—	F01A022610	F01A022612
28	Dark Green	400	—	—	F01A022804	F01A022805	F01A022807	F01A022806	F01A022808	—	F01A022810	F01A022812
30	Grey	400	—	—	F01A023004	F01A023005	F01A023007	F01A023006	F01A023008	—	F01A023010	F01A023012
2-way Female												
12	White	260	—	—	F01A041204	F01A041205	F01A041207	F01A041206	F01A041208	F01A041209	F01A041210	—
14	Green	260	—	—	F01A041404	F01A041405	F01A041407	F01A041406	F01A041408	F01A041409	F01A041410	—
16	Orange	260	—	—	F01A041604	F01A041605	F01A041607	F01A041606	F01A041608	—	F01A041610	F01A041612
18	Red	260	—	—	F01A041804	F01A041805	F01A041807	F01A041806	F01A041808	—	F01A041810	F01A041812
20	Yellow	260	—	—	F01A042004	F01A042005	F01A042007	F01A042006	F01A042008	—	F01A042010	F01A042012
22	Violet	260	—	—	F01A042204	F01A042205	F01A042207	F01A042206	F01A042208	—	F01A042210	F01A042212
24	Blue	260	—	—	F01A042404	F01A042405	F01A042407	F01A042406	F01A042408	—	F01A042410	F01A042412
26	Pink	260	—	—	F01A042604	F01A042605	F01A042607	F01A042606	F01A042608	—	F01A042610	F01A042612
2-way Tiemann												
12	White	400	—	—	F01A081204	F01A081205	F01A081207	F01A081206	F01A081208	F01A081209	F01A081210	—
14	Green	400	—	—	F01A081404	F01A081405	F01A081407	F01A081406	F01A081408	F01A081409	F01A081410	—
16	Orange	400	—	—	F01A081604	F01A081605	F01A081607	F01A081606	F01A081608	—	F01A081610	F01A081612
18	Red	400	—	—	F01A081804	F01A081805	F01A081807	F01A081806	F01A081808	—	F01A081810	F01A081812
20	Yellow	400	—	—	F01A082004	F01A082005	F01A082007	F01A082006	F01A082008	—	F01A082010	F01A082012
22	Violet	400	—	—	F01A082204	F01A082205	F01A082207	F01A082206	F01A082208	—	F01A082210	F01A082212
24	Blue	400	—	—	F01A082404	F01A082405	F01A082407	F01A082406	F01A082408	—	F01A082410	F01A082412
3-way Standard												
16	Orange	400	—	—	F01A101604	F01A101605	F01A101607	F01A101606	F01A101608	—	F01A101610	F01A101612
18	Red	400	—	—	F01A101804	F01A101805	F01A101807	F01A101806	F01A101808	—	F01A101810	F01A101812
20	Yellow	400	—	—	F01A102004	F01A102005	F01A102007	F01A102006	F01A102008	—	F01A102010	F01A102012
22	Violet	400	—	—	F01A102204	F01A102205	F01A102207	F01A102206	F01A102208	—	F01A102210	F01A102212
24	Blue	400	—	—	F01A102404	F01A102405	F01A102407	F01A102406	F01A102408	—	F01A102410	F01A102412
26	Pink	400	—	—	F01A102604	F01A102605	F01A102607	F01A102606	F01A102608	—	F01A102610	F01A102612