



Suction Catheter

■ Features


- Made of DEHP FREE material
- Color-coded connector for easy identification of size
- Soft distal end facilitates comfortable insertion
- Distal end features with side holes or opposite holes
- Available with straight tip & curved tip
- Available with different lengths



■ Position of Aspiration Holes



Standard



Mully



Straight Tip



Curved Tip

■ Straight Tip & Curved Tip

■ Choice of Connectors



Cap-Cone



T-Type



Y-Type



Y-Type Transparent



Plain Type











Funnel Type













Plain Type with Vacuum Interrupter

Standard, Straight Suction Catheter









Size	Color code	Cap-Cone (DEHP FREE)	Cap-Cone (DEHP)	T-Type	Y-Type	Y-Type Transparent	Plain Type	Plain Type with Vacuum Interrupter
6Fr	 Light Green	N02A04061F	N02A040610	N02A02061F	N02A03061F	N02A05061F	N02A01061F	N02A13061F
8Fr	 Light Blue	N02A04081F	N02A040810	N02A02081F	N02A03081F	N02A05081F	N02A01081F	N02A13081F
10Fr	 Black	N02A04101F	N02A041010	N02A02101F	N02A03101F	N02A05101F	N02A01101F	N02A13101F
12Fr	 White	N02A04121F	N02A041210	N02A02121F	N02A03121F	N02A05121F	N02A01121F	N02A13121F
14Fr	 Green	N02A04141F	N02A041410	N02A02141F	N02A03141F	N02A05141F	N02A01141F	N02A13141F
16Fr	 Orange	N02A04161F	N02A041610	N02A02161F	N02A03161F	N02A05161F	N02A01161F	N02A13161F
18Fr	 Red	N02A04181F	N02A041810	N02A02181F	N02A03181F	N02A05181F	N02A01181F	N02A13181F
20Fr	 Yellow	N02A04201F	N02A042010	N02A02201F	N02A03201F	N02A05201F	N02A01201F	N02A13201F

Mully, Straight Suction Catheter

- Improve suction efficiency without damaging the mucous membranes
- Available with different length

Size	Color code	Cap-Cone	Y-Type Transparent	Y-Type Transparent, Soft	Funnel Type	Funnel Type, Soft
4Fr	 Purple	/	/	/	N02A23041F	/
5Fr	 Grey	/	N02A05057F	/	N02A23051F	/
6Fr	 Light Green	N02A04063F	N02A05067F	/	N02A23061F	/
7Fr	 Ivory	/	N02A05077F	/	/	/
8Fr	 Light Blue	N02A04083F	N02A05087F	/	N02A23081F	/
10Fr	 Black	N02A04103F	N02A05107F	N02A26107F	N02A23101F	N02A25101F
12Fr	 White	N02A04123F	N02A05127F	N02A26127F	N02A23121F	N02A25121F
14Fr	 Green	N02A04143F	N02A05147F	N02A26147F	N02A23141F	N02A25141F
16Fr	 Orange	N02A04163F	N02A05167F	N02A26167F	N02A23161F	N02A25161F
18Fr	 Red	N02A04183F	N02A05187F	N02A26187F	N02A23181F	N02A25181F

Mully, Curved Suction Catheter

Size	Color code	Cap-Cone	T-Type	Y-Type	Y-Type Transparent	Plain Type	Plain Type with Vacuum Interrupter	Funnel Type	Funnel Type, Soft
6Fr	 Light Green	N02A08061F	N02A16061F	N02A17061F	N02A19061F	N02A14061F	N02A15061F	/	/
8Fr	 Light Blue	N02A08081F	N02A16081F	N02A17081F	N02A19081F	N02A14081F	N02A15081F	N02A24081F	/
10Fr	 Black	N02A08101F	N02A16101F	N02A17101F	N02A19101F	N02A14101F	N02A15101F	N02A24101F	N02A27101F
12Fr	 White	N02A08121F	N02A16121F	N02A17121F	N02A19121F	N02A14121F	N02A15121F	N02A24121F	N02A27121F
14Fr	 Green	N02A08141F	N02A16141F	N02A17141F	N02A19141F	N02A14141F	N02A15141F	N02A24141F	N02A27141F
16Fr	 Orange	N02A08161F	N02A16161F	N02A17161F	N02A19161F	N02A14161F	N02A15161F	N02A24161F	N02A27161F
18Fr	 Red	N02A08181F	N02A16181F	N02A17181F	N02A19181F	N02A14181F	N02A15181F	N02A24181F	N02A27181F
20Fr	 Yellow	N02A08201F	N02A16201F	N02A17201F	N02A19201F	N02A14201F	N02A15201F	/	/