

Tracheostomy Tube

- Tracheostomy Tube is a hollow tube, with or without a cuff, that is electively inserted directly into the trachea through a surgical incision or with a wire-guided progressive dilatation technique in case of an emergency.
- Lightweight tube made of Medical grade PVC for patient comfort. Full length Radio-opaque line for detection of correct placement. ISO standard connector for universal connection to ventilation equipment. Printed neck plate with size information for easy identification.
- Straps provided in the pack for fixation of the tube. Smooth rounded tip of the Obturator reduces trauma during insertion. Rigid blister pack provides maximum protection for the tube.



■ Cuffed

Size ID (mm)	Cat. No. (DEHP-FREE)	Cat. No. (DEHP)
3.0	A02A05301F	A02A053010
3.5	A02A05351F	A02A053510
4.0	A02A05401F	A02A054010
4.5	A02A05451F	A02A054510
5.0	A02A05501F	A02A055010
5.5	A02A05551F	A02A055510
6.0	A02A05601F	A02A056010
6.5	A02A05651F	A02A056510
7.0	A02A05701F	A02A057010
7.5	A02A05751F	A02A057510
8.0	A02A05801F	A02A058010
8.5	A02A05851F	A02A058510
9.0	A02A05901F	A02A059010
9.5	A02A05951F	A02A059510
10.0	A02A05961F	A02A059610



■ Uncuffed

Size ID (mm)	Cat. No. (DEHP-FREE)	Cat. No. (DEHP)
3.0	A02A06303F	A02A063030
3.5	A02A06353F	A02A063530
4.0	A02A06403F	A02A064030
4.5	A02A06453F	A02A064530
5.0	A02A06503F	A02A065030
5.5	A02A06553F	A02A065530
6.0	A02A06603F	A02A066030
6.5	A02A06653F	A02A066530
7.0	A02A06703F	A02A067030
7.5	A02A06753F	A02A067530
8.0	A02A06803F	A02A068030
8.5	A02A06853F	A02A068530
9.0	A02A06903F	A02A069030
9.5	A02A06953F	A02A069530
10.0	A02A06963F	A02A069630

■ SuctionPlus Tracheostomy Tube

Size ID (mm)	Cat. No. (DEHP-FREE)
6.0	A02C02601F
6.5	A02C02651F
7.0	A02C02701F
7.5	A02C02751F
8.0	A02C02801F
8.5	A02C02851F
9.0	A02C02901F

