



wellead®


Endoscopic Seal

Features:



- **Excellent Self-sealing**
Prevent fluid leakage and reduce risk of cross infection.
- **Good Compatibility**
Suitable for endoscopes with male and female lock fittings.
- **Enhanced Convenience**
A funnel Introducer for easy and safe insertion of instruments up to 3Fr.
Facilitate one-hand manipulation of endoscopes with accessory instruments stabilized by the Seal.


Seal: to mount on the working channel of those endoscopes with male luer lock end.





Introducer: to be inserted into the seal for the passage of guide wire, laser fiber or stone retrieval basket, etc.



Adaptor: available for endoscopes whose working channel has female luer lock end.



Picture	Cat.No.	Contents
	U04D010110	Seal and Introducer
	U04D020110	Seal,Introducer and Adaptor

WELL LEAD MEDICAL CO., LTD.

Address: 47 Guomao Avenue South, Panyu, Guangzhou, China

Tel: +86-20-84758878 Fax: +86-20-84758224 Post Code: 511434

Marketing & Trading Office

Add: RM 805, No 55, 21 Long, Pusan Road, Shanghai, China

Tel: +86-21-58397623 Fax: +86-21-58397615 Post code: 200125

E-mail: info@welllead.com.cn <http://www.welllead.com.cn>

wellead®