


Suction


Yankauer Handle

■ Features

- Made of transparent material for good visualization during suction
- The forming angle of handle is designed to make it more convenient and humanized in clinical use
- The Cross Section is cruciate reinforced rib type, which is designed to make it lightweight and more convenient



With Vent




Without Vent



Standard



Sectional



Crown tip



Round tip



Plain tip

■ Standard Type

Type	Specification	Cat. No.
Round tip	With Vent	N01A020510
Round tip	Without Vent	N01A020610
Plain tip	With Vent	N01A020310
Plain tip	Without Vent	N01A020410

■ Sectional Type

Type	Specification	Cat. No.
Round tip	With Vent	N01A021610
Round tip	Without Vent	N01A021510
Crown tip	With Vent	N01A021110
Crown tip	Without Vent	N01A021210
Plain tip	With Vent	N01A021410
Plain tip	Without Vent	N01A021310



● Suction Connecting Tubing

■ Features

- Made of transparent material for better visualization
- Striated walls of tube provide superior strength and anti-kinking
- Supplied with universal female connector
- Multiple choices of length
- Available with small connector which can connect with suction catheter

Size	Length	Cat. No.
5mm (3/16")	1.8m (6ft)	N01B010110
5mm (3/16")	3m (10ft)	N01B010130
5mm (3/16")	3.6m (12ft)	N01B010136
6mm (1/4")	1.8m (6ft)	N01B010210
6mm (1/4")	3m (10ft)	N01B010230
6mm (1/4")	3.6m (12ft)	N01B010236
7mm (9/32")	1.8m (6ft)	N01B010310
7mm (9/32")	3m (10ft)	N01B010330
7mm (9/32")	3.6m (12ft)	N01B010336

◆ Suction Connecting Tube with Yankauer Handle

■ Features

- Made of transparent material
- Can be combined with different types of yankauer handle and connecting tube
- To suit various requirements



Size	Length	Type	Cat.No.
5mm (3/16")	1.8m (6ft)	Round tip with vent	N01A031510
5mm (3/16")	1.8m (6ft)	Round tip without vent	N01A031610
5mm (3/16")	3m (10ft)	Round tip with vent	N01A031530
5mm (3/16")	3m (10ft)	Round tip without vent	N01A031630
5mm (3/16")	3.6m (12ft)	Round tip with vent	N01A031536
5mm (3/16")	3.6m (12ft)	Round tip without vent	N01A031636
6mm (1/4")	1.8m (6ft)	Round tip with vent	N01A033510
6mm (1/4")	1.8m (6ft)	Round tip without vent	N01A033610
6mm (1/4")	3m (10ft)	Round tip with vent	N01A033530
6mm (1/4")	3m (10ft)	Round tip without vent	N01A033630
6mm (1/4")	3.6m (12ft)	Round tip with vent	N01A033536
6mm (1/4")	3.6m (12ft)	Round tip without vent	N01A033636
7mm (9/32")	1.8m (6ft)	Round tip with vent	N01A036210
7mm (9/32")	1.8m (6ft)	Round tip without vent	N01A036110
7mm (9/32")	3m (10ft)	Round tip with vent	N01A036230
7mm (9/32")	3m (10ft)	Round tip without vent	N01A036130
7mm (9/32")	3.6m (12ft)	Round tip with vent	N01A036236
7mm (9/32")	3.6m (12ft)	Round tip without vent	N01A036136