


Ambulatory Infusion Pump



ISO13485:2016 CE 0123

EASY

Simple operation

Easy parameter
setting

SAFE

Five-year service
life


Complete alarm
system

SMART

Real-time monitoring
system

Multiple infusion modes

Smart delivery



Medsure 3000 series ambulatory infusion pump supports multiple innovative therapies.

The powerful functions make the operation become smart and convenient. The pump is suitable for labor analgesia, postoperative analgesia, tumor chemo, cancerous analgesia and accurate half-life administration.

--Applicable to several departments such as anesthesiology, obstetrics, oncology, pain management, etc.

SPECIFICATIONS

Type:	3000F	<div>DEHP Free</div>
Cassette sizes:	50mL,100mL,150mL,200mL,300mL	
Power supply:	2 AA batteries; at least infusion 100 hours in 2mL/h	
Occlusion level:	High,medium,low ; range: 20kpa ~ 160kpa	
Infusion parameters:		
Total volume	0.1-999ml	
Loading dose	0.1-30ml	
Continuous rate	0.1-120ml/h	
Maximum dose	1-120ml/h	
PCA dose	0.1-9.9ml	
PCA lockout	1-120min	
Intermittent time	1min-48hr59min	
Intermittent dose	1-999ml	
Dose duration	1s-59min59s	
Timed rate	0.1-120ml/h	
Time Settings	24 hours with minimum step 15minutes	
Infusion accuracy	±6%	
Service life	5 years	
Alarms:		
Occlusion alarm; Air alarm; Cassette detached; Near empty; Empty; Low battery; Depleted battery; Key stuck; Reminder alarm; Motor rate error; System error;Data error;		

For details, please refer to the user manual.