

6ml Nebulizer

- Medication volume: 6ml
- Over 80% of particles are $< 5\mu\text{m}$, and the MMAD is $3\mu\text{m}$
- Suitable for short time nebulization



20ml Nebulizer

- Medication volume: 20ml
- Over 80% of particles are $< 5\mu\text{m}$, and the MMAD is $3\mu\text{m}$
- Suitable for situations requiring longer nebulization



Available to be used with aerosol mask and mouthpiece



Aerosol Mask

- Adult Standard
- Adult Elongated
- Pediatric Standard
- Pediatric Elongated



Aerosol Therapy



WELL LEAD MEDICAL CO., LTD.

Address: 47 Guomao Avenue South, Panyu, Guangzhou, China
Tel: +86-20-84758878 Fax: +86-20-84758224 Post Code: 511434
Marketing & Trading Office
Add: RM 805, No 55, 21 Long, Pusan Road, Shanghai, China
Tel: +86-21-58397623 Fax: +86-21-58397615 Post code: 200125
E-mail: info@welllead.com.cn www.welllead.com


wellead

8ml Nebulizer

- Made of transparent polymer material, latex-free, non-toxic and odorless
- Over 90% of particles are < 5µm and the MMAD is 2µm, which can effectively treat lower respiratory tract diseases for the alveoli
- Atomization rate ≥ 3ml/min
- Drive gas flow is about 4 to 8 L/min
- The residual liquid is less than 0.5ml
- Air sourced by compressor


Patented Design

- Leak-proof even when tilted up to 180°
- Greatly improves the atomization rate




8ml Nebulizer - Infant Type

- Atomized particles can effectively reach the alveoli, trachea and bronchi
- Orange color easy to be identified



Available to be used with aerosol mask or mouthpiece



8ml Nebulizer with Aerosol Mask

8ml Nebulizer with Mouthpiece