

Your Complete Enteral Nutrition Portfolio


wellead

VisiSure™

Video Nasogastrintestinal Feeding Tube

Pioneering visual enteral feeding

VisiSure
Video Nasogastrintestinal
Feeding Tube



PU Feeding Tube
Standard Tip & Weighted Tip

Silicone
Feeding Tube

Silicone
Gastrostomy Tube

PVC
Feeding Tube

PVC
Stomach Tube

Silicone
Stomach Tube

WELL LEAD MEDICAL CO., LTD.

Address: 47 Guomao Avenue South, Panyu, Guangzhou, China
Tel: +86-20-84758878 Fax: +86-20-84758224 Post Code: 511434
Marketing & Trading Office
Add: RM 805, No 55, 21 Long, Pusan Road, Shanghai, China
Tel: +86-21-58397623 Fax: +86-21-58397615 Post code: 200125
E-mail: info@wellead.com.cn www.wellead.com

wellead

wl.Brochure.v1.251108



Intended Use

The device is designed for nasal insertion to deliver nutrients, medications, and/or to perform gastrointestinal decompression or irrigation.

Key Features

- Built-in HD Camera with ultra-wide-angle and anti-reflection treatment
- Made of soft TPU material
- Hydrophilic coated tube and stylet
- Coaxial design supports both feeding and decompression
- ENFIT connection design
- Single use

Clinical Value

- Shortens intubation time
- Improves first-attempt success rate
- Enhances patients' nutritional status
- Reduces unplanned extubation (UEX) and adverse events
- Increases nursing satisfaction

ENFIT Connection

Designed for enteral feeding and compliant with ISO 80360-3, enhances patient safety by preventing tubing misconnections

Dual Function in One Tube


Coaxial design for feeding and gastric decompression

Built-in Camera

High definition
Ultra-wide-angle &
anti-reflection

Made of soft TPU and hydrophilic coating, the tube reduces friction and irritation, ensuring safe and comfortable placement

Type	Size	
NI Tube	8Fr / 10Fr / 12Fr / 14Fr / 16Fr / 18Fr / 20Fr	
NGI Tube	10+16Fr / 12+18Fr	



Type-C port